INTUIT SNAPFIT™

TRANSDUCER COVER

REFERENCE GUIDE

INDICATIONS FOR USE / INTENDED USE

Protective cover or sheath placed over diagnostic ultrasound transducer / probe / scanhead instruments. The cover allows use of the transducer in scanning and needle guided procedures for body surface, endocavity, and intra-operative diagnostic ultrasound, while helping to prevent transfer of microorganisms, body fluids, and particulate material to the patient and healthcare worker during reuse of the transducer (both sterile and non-sterile covers). The cover also provides a means for maintenance of a sterile field (sterile covers only). CIVCO Poly Ultrasound Transducer Covers are furnished sterile & non-sterile; single use patient / procedure, disposable.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.


WARNING

- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.
- See your system's user guide for reprocessing transducer between use.
- Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.
- Dispose of single-use components as infectious waste.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- Always place a cover over transducer to protect patients and user from cross-contamination.
- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.
- For proper guide attachment, ensure bands are clear of bracket attachment area. Failure to do so could result in guidance inaccuracy.

NOTE: Cover and bands are not made with natural rubber latex.

APPLYING THE COVER

1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Using proper sterile technique, insert transducer into cover following directional indicators. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Push ring to extend cover over transducer and cord. Bands will deploy as cover is extended.
4. Ensure cover is secured with bands, using additional bands, if necessary.
5. Inspect cover to ensure there are no holes or tears.

6. When procedure is complete, remove ring by twisting slightly to open and pull ring from transducer cord.

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COPYRIGHT © 2017 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. INTUIT SNAPFIT IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
MPS GmbH