REUSABLE ENDOCAVITY NEEDLE GUIDE

REFERENCE GUIDE

For use with
GE Healthcare
IC9
transducers

CE
0120

NEEDLE GUIDE

INDICATIONS FOR USE / INTENDED USE

This device provides a mechanical means for performing needle / instrument guided procedures with the use of the diagnostic ultrasound endocavity transducer.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.
- For illustration purposes only, transducer may be shown without a transducer cover. Always place a cover over transducer to protect patients and user from cross-contamination.
- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- Needle guide is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure needle guide is properly cleaned, sterilized or disinfected before each use.
- Do not use if needle guide is damaged or does not fit properly.


NEEDLE PATH VERIFICATION

- Before using component for first time, perform needle path verification to verify system and biopsy guide relationships, as described in system manual.


NOTE: Use appropriate needle length to reach target area.

USING NEEDLE GUIDANCE


1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Secure cover with enclosed bands.
4. Inspect cover to ensure there are no holes or tears.

5. With clamp open, align guide with locating features on transducer.

6. Close needle guide around transducer and check for secure fit.

CLEANING, DISINFECTING AND STERILIZING NEEDLE GUIDE

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- Dispose of single-use components as infectious waste. Clean and sterilize reusable components after each use.
- See your system's user guide for reprocessing transducer between use.
- This procedure for reprocessing has been validated for effectiveness and compatibility. Device could be damaged or cross-contaminated due to improper reprocessing.
- Inspect needle guide for damage such as cracks or breakage. If damage is evident, discontinue use of needle guide and contact your ultrasound system support representative.

CLEANING THE NEEDLE GUIDE

1. After each use, remove needle guide from transducer, discarding single-use, disposable components. Remove visible contaminants from needle guide surface using a small, soft instrument brush. Keep needle guide from drying out until complete cleaning can be accomplished.
2. Soak needle guide for minimum of five minutes in neutral pH, low foaming enzymatic detergent, such as ENZOL[®] Enzymatic Detergent (Johnson & Johnson). While immersed, use instrument brush to remove trapped contaminants from surfaces and features. If visible contaminants cannot be easily removed, repeat soaking procedure for an additional five minutes. Remove needle guide from cleaning solution and remove any remaining residue with dry wipe. Follow cleaning solution manufacturer's directions for use and recommendations for concentration.

DISINFECTING AND STERILIZING NEEDLE GUIDE

- High-level disinfect needle guide using CIDEX OPA[®] ortho-Phthalaldehyde Solution (Johnson & Johnson) or equivalent .55% ortho-phthalaldehydebased solution. Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedure.
- High-level disinfect or sterilize needle guide using CIDEX[®] Activated Dialdehyde Solution (Johnson & Johnson) or equivalent 2% glutaraldehyde-based solution, CIDEX Plus[®] (Johnson & Johnson) or equivalent 3.4% glutaraldehyde-based solution, or a hydrogen peroxide-based solution. Needle guide may be sterilized by using steam (autoclave) at 132-138°C for 15 to 30 minutes.

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COPYRIGHT © 2017 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
MPS GmbH