



WE'RE ALL ABOUT PEOPLE.

Dornier *Aries*

Orthopedic Shockwave Therapy With Smart Focus Technology



Dornier *MedTech*



ESWT - Non-invasive Treatment

Extracorporeal shockwave therapy (ESWT) is a non-invasive, in-office treatment requiring no anesthesia. Shockwaves are delivered outside of the body to the area of discomfort utilizing a combination of biofeed-back and “Point and Shoot” method, triggering the body’s own repair mechanisms to initiate natural healing.

Shockwave Pioneers

Pioneered by Dornier MedTech, extracorporeal shockwave therapy (ESWT) has been used successfully for a variety of orthopedic indications worldwide for over twenty years.

Ergonomic Applicator

The lightweight, ergonomic design of the applicator allows for comfortable handling during treatments. The integrated control buttons of the applicator permit convenient single hand operation when adjusting the energy level and shockwave frequency.

Advantage

- Wide power range
- Consistent and repeatable shockwaves
- Long EMSE lifetime
- Intuitive operation via touch screen
- Easily exchangeable applicator
- Transportable - Optional transport case or cart



Orthopedic Shockwave Device



Smart Focus

Variable penetration depths with a single applicator are achieved with the special design of the focus channel. Not requiring gel pads, Dornier *Aries* eliminates the need for cumbersome and time consuming exchange during treatment.

Electromagnetic Shockwave Emitter

The Dornier *Aries* utilizes a flat electromagnetic shockwave emitter (EMSE) unique to the orthopedic industry. This technology has demonstrated great success in shockwave applications over past decades due to its consistent and repeatable shockwave characteristics.

Spectrum of Orthopedic Indications

- Lateral epicondylitis
- Plantar fasciitis
- Shoulder and neck pain
- Achillodynia
- Anterior tibia syndrome
- Patellar tendinopathy
- Lower back pain
- Trigger point





WE'RE ALL ABOUT PEOPLE.

Headquarters

Dornier MedTech GmbH
Argelsrieder Feld 7
82234 Wessling
Germany

Europe

Dornier MedTech Europe GmbH
Argelsrieder Feld 7
D-82234 Wessling
Germany
Phone: +49-8153-888 625
Fax: +49-8153-888 444

Dornier MedTech Italia s.r.l.
Via Arrigo Cavaglieri N. 26
00173 Rome
Italy
Phone: +39-06-72 35 289
Fax.: +39-06-72 35 302

Dornier MedTech Espana
Pedro Y Pons 9-11
Edificio Masters
E-08034 Barcelona
Spain
Phone: +34-93-203 9316
Fax: +34-93-205 3814

Dornier MedTech France SARL
Ferme de St Paul
1545 route d'Epagny
F - 74330 Sillingy
France
Phone.: +33-4-50 22 18 94
Fax: +33-4-50 24 24 53

Dornier MedTech Europe GmbH
Moscow Representative Office
Uliza Mytnaya 3
Office 21
119049 Moscow
Russia
Phone: +7-495-230 29 94
Fax.: +7-495-739 51 33

Asia/Pacific

Dornier MedTech Asia Pte. Ltd.
2 Corporation Road
#05-11/12 Corporation Place
Singapore 618494
Phone: +65-6665-5788
Fax: +65-6665 6711

Japan

Dornier MedTech Japan
Meguro Estate Building
3-8-5 Kami-Ohsaki
Shinagawa-ku, Tokyo 141-0021
Japan
Phone: +81-3-3280 3550
Fax: +81-3-3280 3555



CE 1275